

PHCbi Mid Capacity -86°C Freezers

VALID UNTIL:
31 December 2020


Ref. Code:
SANY-MID

For routine storage of specimens at -86°C, we recommend the MDF-U55V-PE Freezer, a 519-litre unit with three shelves and a cascade refrigeration system.

For critical sample storage we recommend the MDF-DU502VX-PE Freezer, a 528-litre unit with the advanced TwinGuard refrigeration system: a dual cooling system that in combination achieves -86°C and if in the rare event one circuit fails, the remaining circuit maintains at -70°C. This superior cooling performance focuses on reliability and safety and is especially suited for mass storage as found in biobanking. With filter-less airflow over the condenser, ensuring no messy filter cleaning required.

- VIP PLUS vacuum insulation panels – large interior space but with smaller footprint
- Insulated inner doors with gasket sealing prevent cold air leakage, removable for easy cleaning
- Comprehensive alarms including high / low temp, power failure, door check, self-diagnostics etc.
- Enhanced safety with dry remote alarm contacts and lockable door
- Quiet in operation
- 3-year warranty on the MDF-U55V-PE and 5-year warranty on the MDF-DU502VX-PE

**20%
OFF!**



Description	Code	Price (ex GST)	Sale Price (ex GST)
Ultra-low Temperature Freezer, 519L with Cascade Refrigeration	SANYMDF-U55V-PE	\$18,976.00	\$14,500.00
Ultra-low Temperature Freezer, 528L with TwinGuard Dual Cooling Refrigeration System	SANYMDF-DU502VX-PE	\$23,893.00	\$18,995.00

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com