

Priorclave QCS Series Autoclaves

VALID UNTIL:
31 December 2020

Ref. Code:
PRIO-AUTO

Order between October and December and receive 10% discount. Contact us now for our recommendations and quotation.

Manufactured in London, the versatile Priorclave QCS series Autoclaves with cylindrical chambers and 2 shelves (front-loading) or baskets (top-loading) have a model suitable for your sterilisation requirements. Constructed to British Standards and meeting AS/NZS standards, the units have grade 316L stainless steel chambers and are inherently safe with thermal and pressure interlocks. Capacities from 100 to 400 litres and a wide range of optional accessories including printers, pre-cycle vacuum and vacuum cooling, vacuum drying, single(standard) 5 or 10 programs, wandering probe, drain condensers, wire and stainless steel baskets, loading cassettes, trolleys and more.

- Choose from: Front Loading: 100/150/200/320 and 400L or Top Loading: 100 and 150L
- Grade 316 stainless steel chamber and door
- Single handed lever action door mechanism
- USB port, fan forced cooling and timed free steaming
- 10-year door and chamber warranty
- Many options and accessories to select from

QCS Series Autoclaves, front-loading or top-loading models to suit your sterilisation requirements.



**10%
OFF!**

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com