

## Whatman™ Uniflo™ Syringe Filters

### Now extended to include 13mm, 25mm and 30mm syringe filters.

Whatman™ Uniflo™ Syringe Filters are disposable filters designed to provide clean filtrate from small volumes up to 100ml. They are available in a variety of membrane choices with a polypropylene overmold housing.

#### Whatman Uniflo syringe filters are available with:

- 13mm, 25mm or 30mm diameters
- 0.2µm or 0.45µm pore sizes
- Sterile or non-sterile options
- Laser etched printing on the filter for easy identification

#### Membranes include:

- Nylon membrane: Aqueous and/or organic samples, hydrophilic
- Polyethersulfone (PES) membrane: Aqueous samples
- Polytetrafluoroethylene (PTFE) membrane: Organic based samples, hydrophobic membrane
- Polyvinylidene difluoride (PVDF) membrane: Aqueous and/or organic based samples
- Low protein binding membrane
- Hydrophilic polytetrafluorethlene (H-PTFE) membrane: Hydrophilic aqueous and/or aggressive organic solvents, hydrophilic, low protein binding
- Chemical compatibility of membranes and housings

**New &  
On Sale!**



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Whatman Uniflo 30mm 0.2um PES w/GF Prefilter 100/PK	WHAT9923-3002	\$370	\$333
Whatman Uniflo 30mm 0.45um PES w/GF Prefilter 100/PK	WHAT9923-3004	\$370	\$333
Whatman Uniflo 30mm 0.2um PES w/GF Prefilter 500/PK	WHAT9924-3002	\$1,464	\$1,318
Whatman Uniflo 30mm 0.45um PES w/GF Prefilter 500/PK	WHAT9924-3004	\$1,464	\$1,318
Whatman Uniflo 30mm 0.2um PVDF w/GF Prefilter 100/PK	WHAT9925-3002	\$370	\$333
Whatman Uniflo 30mm 0.45um PVDF w/GF Prefilter 100/PK	WHAT9925-3004	\$370	\$333
Whatman Uniflo 30mm 0.2um PVDF w/GF Prefilter 500/PK	WHAT9926-3002	\$1,464	\$1,318
Whatman Uniflo 30mm 0.45um PVDF w/GF Prefilter 500/PK	WHAT9926-3004	\$1,464	\$1,318
Whatman Uniflo 30mm 0.2um PTFE w/GF Prefilter 100/PK	WHAT9927-3002	\$370	\$333
Whatman Uniflo 30mm 0.45um PTFE w/GF Prefilter 100/PK	WHAT9927-3004	\$370	\$333
Whatman Uniflo 30mm 0.2um PTFE w/GF Prefilter 500/PK	WHAT9928-3002	\$1,464	\$1,318
Whatman Uniflo 30mm 0.45um PTFE w/GF Prefilter 500/PK	WHAT9928-3004	\$1,464	\$1,318
Whatman Uniflo 30mm 0.2um Nylon w/GF Prefilter 100/PK	WHAT9929-3002	\$370	\$333
Whatman Uniflo 30mm 0.45um Nylon w/GF Prefilter 100/PK	WHAT9929-3004	\$370	\$333
Whatman Uniflo 30mm 0.2um Nylon w/GF Prefilter 500/PK	WHAT9930-3002	\$1,464	\$1,318
Whatman Uniflo 30mm 0.45um Nylon w/GF Prefilter 500/PK	WHAT9930-3004	\$1,464	\$1,318
Whatman Uniflo 25mm 0.45um PVDF 100/PK	WHAT9917-2504	\$246	\$221
Whatman Uniflo 25mm 0.45um Nylon 100/PK	WHAT9918-2504	\$190	\$171
Whatman Uniflo H-PTFE 25/0.2um 100/PK	WHAT9920-2502	\$188	\$169
Whatman Uniflo H-PTFE 25/0.2um 500/PK	WHAT9921-2502	\$772	\$695
Whatman Uniflo H-PTFE 25/0.45um 100/PK	WHAT9920-2504	\$188	\$169
Whatman Uniflo H-PTFE 25/0.45um 500/PK	WHAT9921-2504	\$772	\$695
Whatman Uniflo H-PTFE 13/0.2um 100/PK	WHAT9920-1302	\$175	\$158
Whatman Uniflo H-PTFE 13/0.2um 500/PK	WHAT9921-1302	\$667	\$601
Whatman Uniflo H-PTFE 13/0.45um 100/PK	WHAT9920-1304	\$175	\$158
Whatman Uniflo H-PTFE 13/0.45um 500/PK	WHAT9921-1304	\$667	\$601

## New Whatman Uniflo Syringe Filters on Sale!

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Whatman Uniflo 13mm 0 22 PVDF 500/PK	WHAT9909-1302	\$947	\$852
Whatman Uniflo 13mm 0 45 PVDF 500/PK	WHAT9909-1304	\$947	\$852
Whatman Uniflo 13mm 0 22 Nylon 500/PK	WHAT9910-1302	\$927	\$834
Whatman Uniflo 13mm 0 45 Nylon 500/PK	WHAT9910-1304	\$927	\$834
Whatman Uniflo 13mm 0 22 PTFE 500/PK	WHAT9911-1302	\$947	\$852
Whatman Uniflo 13mm 0 45 PTFE 500/PK	WHAT9911-1304	\$947	\$852
Whatman Uniflo 13mm 0 45 PES 500/PK	WHAT9912-1304	\$927	\$834
Whatman Uniflo 25mm 0 22 PVDF 500/PK	WHAT9909-2502	\$1,205	\$1,084
Whatman Uniflo 25mm 0 45 PVDF 500/PK	WHAT9909-2504	\$1,185	\$1,066
Whatman Uniflo 25mm 0 22 Nylon 500/PK	WHAT9910-2502	\$1,185	\$1,066
Whatman Uniflo 25mm 0 45 Nylon 500/PK	WHAT9910-2504	\$1,185	\$1,066
Whatman Uniflo 25mm 0 22 PES 500/PK	WHAT9912-2502	\$1,185	\$1,066
Whatman Uniflo 25mm 0 45 PES 500/PK	WHAT9912-2504	\$1,185	\$1,066
Whatman Uniflo 13mm 0 22 PES 500/PK	WHAT9912-1302	\$927	\$834
Whatman Uniflo 25mm 0 22 PTFE 500/PK	WHAT9911-2502	\$1,228	\$1,105
Whatman Uniflo 25mm 0 45 PTFE 500/PK	WHAT9911-2504	\$1,228	\$1,105
Whatman Uniflo 25mm 0 22 PVDF S 45/PK	WHAT9913-2502	\$96	\$86
Whatman Uniflo 25mm 0 45 PVDF S 45/PK	WHAT9913-2504	\$96	\$86
Whatman Uniflo 25mm 0 2 PES S 45/PK	WHAT9914-2502	\$96	\$86
Whatman Uniflo 25mm 0 45 PES S 45/PK	WHAT9914-2504	\$96	\$86
Whatman Uniflo 25mm 0 2 PES S 200/PK	WHAT9915-2502	\$426	\$383
Whatman Uniflo 13mm 0 2 PES S 100/PK	WHAT9916-1302	\$213	\$192
Whatman Uniflo 13mm 0 45 PES S 100/PK	WHAT9916-1304	\$213	\$192
Whatman Uniflo 25mm 0 45 PES S 200/PK	WHAT9915-2504	\$443	\$398

# Whatman Uniflo Syringe Filters

Consistent performance. Maximum value.



Contact us for more information  
and to place your orders

**Bio-Strategy Ltd.**

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

[shop.bio-strategy.com](http://shop.bio-strategy.com)



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 June 2021

Ref. Code:  
WHAT-UNIFLO

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.