

## Axygen Microcentrifuge Tubes

### Choice of 0.6ml, 1.5ml, 1.7ml or 2ml MaxyClear, coloured, Maxymum Recovery® or standard surface tubes.

Axygen microcentrifuge tubes are available in MaxyClear or an assortment of colours including an amber option. Produced from virgin polypropylene with either Axygen's Maxymum Recovery® (L-C) for minimal sample retention or the standard surface. Snaplock microcentrifuge tubes have frosted cap surfaces and frosted side panels for easy labelling. Each cap centre features a thin section piercing port for easy syringe sampling.

- MaxyClear Homopolymer 0.6ml, 1.5ml and 2ml microcentrifuge tubes are boil-proof
- Co-polymer 1.7ml tubes offer exceptional strength and chemical resistance
- Most tubes are autoclavable
- Maximum G-force is 14,000 RCF
- RNase-/DNase-free and nonpyrogenic
- Sterile option also available

**Special New  
Year Prices**



Description	Pack Size	Code	Price (excl. GST)	Sale Price (excl. GST)
Microtube 0.6ml MaxyClear Snaplock MR	500 Items	AXYGMCT-060-L-C	<b>\$33</b>	<b>\$25</b>
Microtube 0.6ml MaxyClear Snaplock	1000 Items	AXYGMCT-060-C	<b>\$43</b>	<b>\$31</b>
Microtube 0.6ml Violet Snaplock	1000 Items	AXYGMCT-060-V_U	<b>\$45</b>	<b>\$31</b>
Microtube 0.6ml Assorted Colour Snaplock	1000 Items	AXYGMCT-060-A	<b>\$47</b>	<b>\$36</b>
Microtube 0.6ml Blue Snaplock	1000 Items	AXYGMCT-060-B_U	<b>\$45</b>	<b>\$31</b>
Microtube 1.5ml MaxyClear Snaplock	500 Items	AXYGMCT-150-C	<b>\$21</b>	<b>\$16</b>
Microtube 1.5ml MaxyClear Snaplock MR	250 Items	AXYGMCT-150-L-C	<b>\$18</b>	<b>\$15</b>
Microtube 1.7ml MaxyClear Snaplock	500 Items	AXYGMCT-175-C	<b>\$19</b>	<b>\$14</b>
Microtube 1.7ml MaxyClear Snaplock MR	250 Items	AXYGMCT-175-L-C	<b>\$18</b>	<b>\$15</b>
Microtube 1.7ml Amber Snaplock	500 Items	AXYGMCT-175-X	<b>\$26</b>	<b>\$21</b>
Microtube 1.7ml Red Snaplock	500 Items	AXYGMCT-175-R	<b>\$26</b>	<b>\$18</b>
Microtube 1.7ml Blue Snaplock	500 Items	AXYGMCT-175-B	<b>\$26</b>	<b>\$18</b>
Microtube 1.7ml Assorted Colour Snaplock	500 Items	AXYGMCT-175-A	<b>\$27</b>	<b>\$21</b>
Microtube 2ml MaxyClear Snaplock	500 Items	AXYGMCT-200-C	<b>\$23</b>	<b>\$17</b>
Microtube 2ml MaxyClear No Cap	5000 Items	AXYGMCT-200-NC	<b>\$232</b>	<b>\$160</b>
Microtube 2ml MaxyClear Snaplock MR	2500 Items	AXYGMCT-200-L-C	<b>\$177</b>	<b>\$154</b>

Contact us for more information  
and to place your orders

#### Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 March 2021

Ref. Code:  
AXYG-SNAP

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.