

## Stock Up and Save on Azlon's HDPE and PP Storage Bottles!

### Ring in the New Year with these super specials.

#### Stock up on your bottles now and save!

A range of HDPE and PP bottles, useful in any laboratory.

- Ideal for a wide range of solid or liquid storage, media preparation or sampling applications
- Range of volumes in HDPE or PP, round or square, wide or narrow neck to suit your needs
- Suitable for food contact applications in accordance with EU and U.S. FDA regulations

**SAVE  
TODAY!**



Description	Pack Size	Code	Price (excl. GST)	Sale Price (excl. GST)
Bottle 25ml HDPE Wide Mouth Square Plug Seal Blue Screw Cap	10 Items	AZLOBGE212P	\$21	<b>\$14</b>
Bottle 250ml HDPE Wide Mouth Square Plug Seal Blue Screw Cap	10 Items	AZLOBGE230P	\$31	<b>\$22</b>
Bottle 100ml HDPE Wide Neck PP Screw Cap	10 Items	AZLOBLH0100P	\$10	<b>\$7</b>
Bottle 125ml HDPE Narrow Neck PP Screw Cap	10 Items	AZLOBNH0125PN	\$8	<b>\$6</b>
Bottle 150ml HDPE Narrow Neck PP Screw Cap	10 Items	AZLOBNH0150PN	\$9	<b>\$6</b>
Bottle 250ml HDPE Narrow Neck PP Screw Cap	10 Items	AZLOBNH0250P	\$15	<b>\$10</b>
Bottle 500ml HDPE Narrow Neck PP Screw Cap	10 Items	AZLOBNH0500P	\$19	<b>\$13</b>
Bottle 250ml HDPE Wide Mouth PP Screw Cap	10 Items	AZLOBWH0250PN	\$17	<b>\$12</b>
Bottle 500ml HDPE Wide Mouth PP Screw Cap	10 Items	AZLOBWH0500PN	\$25	<b>\$17</b>
Bottle 1000ml HDPE Wide Mouth PP Screw Cap	5 Items	AZLOBWH1000PN	\$27	<b>\$19</b>
Bottle 2000ml HDPE Wide Mouth PP Screw Cap	5 Items	AZLOBWH2000P	\$33	<b>\$23</b>
Bottle 250ml PP Wide Mouth Heavy Duty Screw Cap	10 Items	AZLOBWP0250P	\$21	<b>\$15</b>
Bottle 250ml PP Wide Mouth Graduated Screw Cap	10 Items	AZLOBWP0250PGRAD	\$50	<b>\$36</b>
Bottle 1000ml PP Wide Mouth Heavy Duty Screw Cap	5 Items	AZLOBWP1000P	\$27	<b>\$19</b>

Contact us for more information  
and to place your orders

#### Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 March 2021

Ref. Code:  
AZLO-HDPE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

# Azlon®

## AMBER POLYPROPYLENE STORAGE BOTTLES

Now meet US pharmacopoeia 39 requirements  
for light transmission

**Amber polypropylene storage bottles** are designed to minimise UV light transmissions and to protect light-sensitive substances.

They now conform to USP 39 requirements so you can be sure that your light-sensitive substances are safe.



CONTACT SCILABWARE TO DISCUSS  
OPTIONS AND PRICING

Visit our website for more information  
[www.scilabware.com](http://www.scilabware.com)

Email us with your questions  
[azlon@scilabware.com](mailto:azlon@scilabware.com)

**SciLabware Limited** Unit 4, Riverside 2,  
Campbell Road, Stoke on Trent, ST4 4RJ, UK.  
Tel: +44 (0)1782 444406 Fax: +44 (0)1782 940436



# SciLabware

# Azlon®

## STORAGE BOTTLES IN BULK QUANTITIES

**Did you know... Azlon® storage bottles are available in bulk quantities and uncapped for a swifter packaging process...**

**Azlon® storage bottles are available in many shapes and sizes, with a leakproof, linerless design to ensure a secure closure.**

Manufactured using the highest quality resins, and with a range of various polymer options, you'll find the perfect solution whatever your application.

### **Leakproof**

High quality moulded threads and seals ensure secure closure.

### **Linerless**

Eliminate contamination from loose cap insert liners.

### **Base mouldings**

Include recycling code, bottle volume and registration ramp.



# SciLabware

**CONTACT SCILABWARE TO DISCUSS  
OPTIONS AND PRICING**

Visit our website for more information  
[www.scilabware.com](http://www.scilabware.com)

Email us with your questions  
[azlon@scilabware.com](mailto:azlon@scilabware.com)

**SciLabware Limited** Unit 4, Riverside 2,  
Campbell Road, Stoke on Trent, ST4 4RJ, UK.  
Tel: +44 (0)1782 444406 Fax: +44 (0)1782 940436



# STORAGE BOTTLES IN BULK QUANTITIES

Product Code	Bottle Polymer	Cap Polymer	Neck Type	Capacity (ml)	Height (mm)	O.D. (mm)	Neck I.D (mm)	Case Quantity	Half Pallet Quantity	Full Pallet Quantity
BNP0030BULK	PP	PP	NARROW	30	66	32	13	1200	12000	24000
BNP0060BULK	PP	PP	NARROW	60	80	40	13	800	8000	16000
BNP0100BULK	PP	PP	NARROW	100	89	50	18	420	4200	8400
BNP0125BULK	PP	PP	NARROW	125	102	50	18	350	3500	7000
BNP0150BULK	PP	PP	NARROW	150	116	50	18	330	3300	6600
BNP0250BULK	PP	PP	NARROW	250	129	60	18	200	2000	4000
BNP0500BULK	PP	PP	NARROW	500	167	75	21	110	1100	2200
BNP1000BULK	PP	PP	NARROW	1000	213	90	27	60	600	1200
BNP2000BULK	PP	PP	NARROW	2000	240	120	27	30	300	600
BNH0030BULK	HDPE	PP	NARROW	30	66	32	13	1200	12000	24000
BNH0060BULK	HDPE	PP	NARROW	60	80	40	13	800	8000	16000
BNH0100BULK	HDPE	PP	NARROW	100	89	50	18	420	4200	8400
BNH0125BULK	HDPE	PP	NARROW	125	102	50	18	350	3500	7000
BNH0150BULK	HDPE	PP	NARROW	150	116	50	18	330	3300	6600
BNH0250BULK	HDPE	PP	NARROW	250	129	60	18	200	2000	4000
BNH0500BULK	HDPE	PP	NARROW	500	167	75	21	110	1100	2200
BNH1000BULK	HDPE	PP	NARROW	1000	213	90	27	60	600	1200
BNH2000BULK	HDPE	PP	NARROW	2000	240	120	27	30	300	600
BNL0030BULK	LDPE	PP	NARROW	30	66	33	13	1200	12000	24000
BNL0060BULK	LDPE	PP	NARROW	60	80	40	13	800	8000	16000
BNL0100BULK	LDPE	PP	NARROW	100	89	50	18	420	4200	8400
BNL0125BULK	LDPE	PP	NARROW	125	102	50	18	350	3500	7000
BNL0150BULK	LDPE	PP	NARROW	150	116	50	18	330	3300	6600
BNL0250BULK	LDPE	PP	NARROW	250	129	62	18	200	2000	4000
BNL0500BULK	LDPE	PP	NARROW	500	167	73	21	110	1100	2200
BNL1000BULK	LDPE	PP	NARROW	1000	213	92	27	60	600	1200
BNL2000BULK	LDPE	PP	NARROW	2000	240	121	27	30	300	600
BWP0030BULK	PP	PP	WIDE	30	70	32	22	1200	12000	24000
BWP0060BULK	PP	PP	WIDE	60	82	40	22	800	8000	16000
BLP0100BULK	PP	PP	WIDE	100	89	51	32	420	4200	8400
BLP0125BULK	PP	PP	WIDE	125	104	51	32	350	3500	7000
BLP0150BULK	PP	PP	WIDE	150	116	51	32	330	3300	6600
BWP0250BULK	PP	PP	WIDE	250	140	60	37	200	2000	4000
BWP0500BULK	PP	PP	WIDE	500	161	75	46	110	1100	2200
BWP1000BULK	PP	PP	WIDE	1000	205	90	55	60	600	1200
BWP2000BULK	PP	PP	WIDE	2000	242	120	55	30	300	600
BWH0030BULK	HDPE	PP	WIDE	30	70	32	22	1200	12000	24000
BWH0060BULK	HDPE	PP	WIDE	60	82	40	22	800	8000	16000
BLH0100BULK	HDPE	PP	WIDE	100	89	51	32	420	4200	8400
BLH0125BULK	HDPE	PP	WIDE	125	104	51	32	350	3500	7000
BLH0150BULK	HDPE	PP	WIDE	150	116	51	32	330	3300	6600
BWH0250BULK	HDPE	PP	WIDE	250	140	60	37	200	2000	4000
BWH0500BULK	HDPE	PP	WIDE	500	161	75	46	110	1100	2200
BWH1000BULK	HDPE	PP	WIDE	1000	205	90	55	60	600	1200
BWH2000BULK	HDPE	PP	WIDE	2000	242	120	55	30	300	600