

## Esco Airstream E Class II Microbiological Safety Cabinet, 1.2m

**Free adjustable height stand, UV light and 2 power sockets (valued at \$1,260).**

**Ensuring your safety in the laboratory with operator, product and environment protection.**

With over forty years manufacturing experience in laminar flow, the AC2-4E8 is designed with your health and safety in mind. The ergonomic sloping front combined with tempered glass side panels allows entry of more natural light. The manual sash provides trouble-free smooth operation with automatic activation of lighting and airflow when set to the correct working height.

- Interlocked UV light with timer – no accidental burns
- Airflow monitoring and alarms – cabinet operating status clearly displayed
- Multi-piece stainless steel work surface – easy to remove and clean
- 53.5dB(A): low noise level for a more comfortable working environment
- 3-year warranty

**15%  
OFF!**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Airstream E Class II Type A2 Microbiological Safety Cabinet, 1.2m	ESCOAC2-4E8	\$13,031	<b>\$10,755</b>



Contact us for more information  
and to place your orders

**Bio-Strategy Ltd.**

T: 0800 34 24 66  
E: sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 March 2021

Ref. Code:  
ESCO-SAFE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.