

Heidolph Hei-Torque Core Overhead Stirrer

Constant speed even under changing loads!

Price valid until stock lasts!

- Mix up to 25L
- Speed range up to 2,000rpm
- Viscosity up to 10,000mPas
- Control knob for quick settings
- Timer
- No tools required to change an impeller
- Germany quality: 3-year warranty!


Do you need to mix larger volumes?
Do you process high viscous fluids?

Hei-TORQUE Value: The reliable overhead stirrer for standard applications. Clearly laid out display and great ease of operation. They perform stirring tasks quickly and reliably.

Hei-TORQUE Precision models are ideal for demanding tasks that have to be reproducible and documentable.

- Forceful stirring in three performance classes up to 400Ncm (for up to 250,000mPas)
- Mix up to 100L
- Speed range up to 2,000rpm even mixing liquids with high viscosity

**30%
OFF!**



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Hei-Torque Core	HEID501-60410-00	\$1,563	\$1,313
BR 13 Square-Blade Impeller, for mixing tasks with little or average viscosity	HEID509-13000-00	\$114	\$79
VISCO JET® 80 mm, (Stainless Steel), the only impeller world wide capable of completely mixing larger quantities of high-viscosity liquids and gels	HEID509-16080-00	\$1,208	\$955
PR 30 Pitched-Blade Impeller (Stainless Steel), for mixing tasks with medium or high viscosity; for applications which require average or high speed, excellent mixing properties for homogenisation and suspensions	HEID509-30000-00	\$127	\$103
PR 39 Pitched-Blade Impeller(PTFE), for mixing tasks with medium or high viscosity; for applications which require average or high speed. Excellent mixing properties for homogenization and suspensions	HEID509-39000-10	\$485	\$414

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2021

Ref. Code:
HEID-CORE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.