

Miele Professional PG8593 Washer Disinfector with DryPlus HEPA Filtered Air Drying

**Free external DOS K85 dispensing/
dosing pump for detergent or
neutraliser (valued at \$1,189).**

Laboratory underbench/free-standing washer with stainless steel interior combining powerful pump with integrated heating elements enabling safe washing cycles up to 93°C. The unit also features spray arm and spray pressure monitoring with alerts to advise you if a blockage or pressure deviation has occurred.

- 15 programs including new "Mini" program for lightly soiled glassware
- Internal dispensing pump for liquid detergent or liquid neutraliser
- Compact: 600mm x 600mm x 850mm (820mm with lid removed)
- Many other options and accessories to select from

Description	Code	Price (excl. GST)
Miele Professional PG8593 Washer Disinfector with DryPlus HEPA Filtered Air Drying	MIELG8593-Package1	\$19,583



Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2021

Ref. Code:
MIEL-PG8593

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.