

## PHCbi MDF-DU702VX-PE: 729 Litre Upright Freezer

### One of the safest ULT freezers on the market today!

The best redundancy you can have! If one refrigeration circuit fails, the other circuit will maintain -70°C.

- **TwinGuard® parallel compressors:** for increased product storage safety
- **VIP Plus insulation:** provides typically 24% more storage capacity compared to conventionally insulated freezers with the same footprint
- **User friendly colour graphical control panel**
- **Intelligent ECO mode:** optimising compressor operation to match the freezer's running condition and safe electricity costs
- **Filter-less condenser design:** eliminates the need for routine filter cleaning

**Save Over  
\$6,000!**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
PHCbi 729 Litre Upright ULT freezer with TwinGuard Refrigeration	SANY MDF-DU702VX-PE	<b>\$26,154</b>	<b>\$19,775</b>



Contact us for more information  
and to place your orders

#### Bio-Strategy Ltd.

T: 0800 34 24 66  
E: sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 March 2021

Ref. Code:  
SANY-MDF

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.