

R-Biopharm RIDA®CUBE SCAN

The handy analyser for enzymatic single-testing.

- Allows single sample testing in small wine, beverage and food testing facilities
- The photometric system allows biochemistry testing, covering all enzymatic and colorimetric assays for the detection of organic acids (e.g. lactic acid), sugars (e.g. glucose) or other food components (e.g. sulfite)
- Single-test cartridges with ready-to-use reagents
- RFID-card containing all test settings
- Usually only one pipetting step for the sample, results in 15 minutes
- Already have a spectrophotometer - we can supply our Roche "Yellow Line" or Enzytec liquid stable enzymatic kits to run on those – special offer 10% discount, please enquire for details

Please speak to your local Bio-Strategy Sales Specialist for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

**10%
OFF!**



Description	Pack Size	Code	Price (excl. GST)	Sale Price (excl. GST)
Rida®Cube Scan 340 / 580nm Instrument	1 Item	RBIOZRCS0580	\$4,000	\$3,600
Rida®Cube Scan Lactose/D-Galactose Enzymatic Test 340nm	32 Test Kit	RBIORCS4110	\$257	\$231
Rida®Cube Scan D-Galactose Enzymatic Test 340nm	32 Test Kit	RBIORCS4120	\$228	\$205
Rida®Cube Scan D-Glucose Enzymatic Test 340nm	32 Test Kit	RBIORCS4140	\$199	\$179
Rida®Cube Scan D-Glucose/D-Fructose Enzymatic Test 340nm	32 Test Kit	RBIORCS4160	\$238	\$214
Rida®Cube Scan Sucrose/D-Glucose Enzymatic Test 340nm	32 Test Kit	RBIORCS4180	\$241	\$217
Rida®Cube Scan Sucrose/D-Glucose/D-Fructose Enzymatic Test 340nm	32 Test Kit	RBIORCS4190	\$309	\$278
Rida®Cube Scan Acetic Acid Enzymatic Test 340nm	32 Test Kit	RBIORCS4226	\$315	\$284
Rida®Cube Scan D/L-Lactic Acid Enzymatic Test 340nm	32 Test Kit	RBIORCS4240	\$363	\$327
Rida®Cube Scan L-Lactic Acid Enzymatic Test 340nm	32 Test Kit	RBIORCS4260	\$257	\$231
Rida®Cube Scan L-Malic Acid Enzymatic Test 340nm	32 Test Kit	RBIORCS4280	\$225	\$203
Rida®Cube Scan Ethanol Enzymatic Test 340nm	32 Test Kit	RBIORCS4340	\$268	\$241
Rida®Cube Scan SO ₂ -Total (Total Sulfite) Enzymatic Test 340nm	32 Test Kit	RBIORCS4600	\$206	\$185
Rida®Cube Scan Ammonia Enzymatic Test 340nm	32 Test Kit	RBIORCS4390	\$337	\$303
Rida®Cube Scan SO ₂ -Free (Free Sulfite) Enzymatic Test 340nm	32 Test Kit	RBIORCS4610	\$274	\$246

Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2021

Ref. Code:
RBIO-CUBE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.