


We back you up — Consumables for Preparation of COVID-19 Tests



Your source of high quality plates and sealing films for RT-PCR machines

Brand	Model	Compatible plates	Product Description	Case Qty
Applied Biosystems	QuantStudio™ qPCR Systems	PCR-96-AB-C	200 µl, 96-well, elevated skirt	10/pk, 5pk/cs
		PCR-96M2-HS-C	200 µl, 96-well, half skirt, single notch corner	10/pk, 5pk/cs
	ViiA 7™ Real-Time PCR System	PCR-96-LP-AB-C (fast reaction)	100 µl, 96-well PCR microplate, low profile	5 sleeves/pk, 5pk/cs
		PCR-384M2-C	55 µl, 384-well, skirted	10/pk, 5pk/cs
Roche	LightCycler® 480 Real-Time PCR System	PCR-96-LC480-W-NF	100 µl, 96-well PCR microplate, white, no sealing film	10/pk, 5pk/cs
		PCR-384-LC480-W-NF	100 µl, 96-well PCR microplate, white, no sealing film	10/pk, 5pk/cs
All brands	All models	UC-500	70 µm Ultra Clear Pressure Sensitive Sealing Film for Real Time PCR, Nonsterile	100/pk, 5pk/cs



Learn more about
Axygen® PCR Microplates



Learn more about
Axygen® sealing films

Your Source of High Quality Tips for Automated Liquid Handlers




Axygen® automation tips are high quality alternatives to more expensive tips sold by automation equipment manufacturers. We specialize in manufacturing tips for all major automation platforms, providing equal—or superior—performance to the tips supplied by automation equipment manufacturers.

[Learn more](#)

Your source of high quality deep well plate and reservoirs for Automation solutions

Axygen® assay and deep well plates are used for sample collection, *in vitro* growth chamber, and long-term storage. Made of polypropylene. They offer excellent chemical resistance and will withstand temperatures down to -80°C.



▶ **Axygen® Deep Well**

[Learn more](#)

in **New Zealand** please contact
Bio-Strategy Ltd

CORNING

FALCON

AXYGEN

PYREX

T: 0800 34 24 66 | E: sales.nz@bio-strategy.com
www.bio-strategy.com | shop.bio-strategy.com